

FGD3 antibody

Product Information

Catalog No.:	FNab03085
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Promotes the formation of filopodia. May activate CDC42, a member of the Ras-like family of Rho-and Rac proteins, by exchanging bound GDP for free GTP. Plays a role in regulating the actin cytoskeleton and cell shape(By similarity).

Immunogen information


Immunogen:	FYVE, RhoGEF and PH domain containing 3
Synonyms:	FYVE, RhoGEF and PH domain-containing protein 3 Zinc finger FYVE domain-containing protein 5 FGD3 ZFYVE5
Observed MW:	79 kDa
Uniprot ID :	Q5JSP0

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab03085 (FGD3 antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab03085 (FGD3 antibody) at dilution of 1:500