

## FGF12 antibody

### Product Information

Catalog No.:	FNab03088
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in nervous system development and function. Promote neuronal excitability by elevating the voltage dependence of neuronal sodium channel SCN8A fast inactivation.

### Immunogen information


Immunogen:	fibroblast growth factor 12
Synonyms:	Fibroblast growth factor 12 (FGF-12) Fibroblast growth factor homologous factor 1 (FHF-1) Myocyte-activating factor FGF12 FGF12B FHF1
Observed MW:	20-25 kDa
Uniprot ID :	P61328

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab03088(FGF12 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03088(FGF12 antibody) at dilution of 1:1000