

## FGF17 antibody

### Product Information

Catalog No.:	FNab03092
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Plays an important role in the regulation of embryonic development and as signaling molecule in the induction and patterning of the embryonic brain. Required for normal brain development.

### Immunogen information

Immunogen:	fibroblast growth factor 17
Synonyms:	Fibroblast growth factor 17 (FGF-17) FGF17
Observed MW:	33 kDa
Uniprot ID :	O60258

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab03092( FGF17 Antibody) at dilution of 1:600