

FGFBP1 antibody

Product Information

Catalog No.:	FNab03096
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a carrier protein that release fibroblast-binding factors(FGFs) from the extracellular matrix(EM) storage and thus enhance the mitogenic activity of FGFs. Enhances FGF2 signaling during tissue repair, angiogenesis and in tumor growth.

Immunogen information

Immunogen:	fibroblast growth factor binding protein 1
Synonyms:	Fibroblast growth factor-binding protein 1 (FGF-BP, FGF-BP1, FGF-binding protein 1, FGFBP-1) 17 kDa heparin-binding growth factor-binding protein (17 kDa HBGF-binding protein, HBp17) FGFBP1 FGFBP HBP17
Observed MW:	26 kDa
Uniprot ID :	Q14512

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03096(FGFBP1 antibody) at dilution of 1:50