

## FGFBP3 antibody

### Product Information

Catalog No.:	FNab03098
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Heparin-binding protein which binds to FGF2, prevents binding of FGF2 to heparin and probably inhibits immobilization of FGF2 on extracellular matrix glycosaminoglycans, allowing its release and subsequent activation of FGFR signaling which leads to increased vascular permeability.

### Immunogen information


Immunogen:	fibroblast growth factor binding protein 3
Synonyms:	Fibroblast growth factor-binding protein 3 (FGF-BP3, FGF-binding protein 3, FGFBP-3) FGFBP3 C10orf13
Observed MW:	60-65 kDa
Uniprot ID :	Q8TAT2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon tissue slide using FNab03098(FGFBP3 Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab03098(FGFBP3 Antibody) at dilution of 1:300