

FHL5 antibody

Product Information

Catalog No.:	FNab03113
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in the regulation of spermatogenesis. Stimulates CREM transcriptional activity in a phosphorylation-independent manner.

Immunogen information

Immunogen:	four and a half LIM domains 5
Synonyms:	Four and a half LIM domains protein 5 (FHL-5) Activator of cAMP-responsive element modulator in testis (Activator of CREM in testis) FHL5 ACT
Observed MW:	33 kDa
Uniprot ID :	Q5TD97

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab03113(FHL5 Antibody) at dilution of 1:2000