

## FBN2-Specific antibody

### Product Information

Catalog No.:	FNab03116
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	fibrillin 2
Synonyms:	Fibrillin-2 Placensin FBN2
Observed MW:	290 kDa
Uniprot ID :	P35556

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab03116(FBN2-Specific Antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03116(FBN2-Specific antibody) at dilution of 1:500