

FBLN5 antibody

Product Information


Catalog No.:	FNab03126
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8H8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	fibulin 5
Synonyms:	Fibulin-5 (FIBL-5) Developmental arteries and neural crest EGF-like protein (Dance) Urine p50 protein (UP50) FBLN5 DANCE
Observed MW:	64-66 kDa
Uniprot ID :	Q9UBX5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab03126(FBLN5 antibody) at dilution of 1:50



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using FNab03126 (FBLN5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab03126(FBLN5 antibody) at dilution of 1:500