

FBLN1 antibody

Product Information


Catalog No.:	FNab03127
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	fibulin 1
Synonyms:	Fibulin-1 (FIBL-1) FBLN1
Observed MW:	77 kDa
Uniprot ID :	P23142

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab03127(Fibulin-1 antibody at dilution of 1:1000