

FKBP3 antibody

Product Information

Catalog No.:	FNab03140
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


FK506-and rapamycin-binding proteins(FKBPs) constitute a family of receptors for the two immunosuppressants which inhibit T-cell proliferation by arresting two distinct cytoplasmic signal transmission pathways. PPIases accelerate the folding of proteins.

Immunogen information

Immunogen:	FK506 binding protein 3, 25 kDa
Synonyms:	Peptidyl-prolyl cis-trans isomerase FKBP3 (PPIase FKBP3) 25 kDa FK506-binding protein (25 kDa FKBP, FKBP-25) FK506-binding protein 3 (FKBP-3) Immunophilin FKBP25 Rapamycin-selective 25 kDa immunophilin Rotamase FKBP3 FKBP25
Observed MW:	25 kDa
Uniprot ID :	Q00688

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab03140(FKBP3 antibody) at dilution of 1:1000