

FKBP7 antibody

Product Information

Catalog No.:	FNab03148
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


PPIases accelerate the folding of proteins during protein synthesis.

Immunogen information


Immunogen:	FK506 binding protein 7
Synonyms:	Peptidyl-prolyl cis-trans isomerase FKBP7 (PPIase FKBP7) 23 kDa FK506-binding protein (23 kDa FKBP, FKBP-23) FK506-binding protein 7 (FKBP-7) Rotamase FKBP7 FKBP23
Observed MW:	30 kDa
Uniprot ID :	Q9Y680

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab03148(FKBP7 antibody) at dilution of 1:50



human lung tissue were subjected to SDS PAGE followed by western blot with FNab03148(FKBP7 antibody) at dilution of 1:500