

FKBP9 antibody

Product Information

Catalog No.:	FNab03150
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


PPIases accelerate the folding of proteins during protein synthesis.

Immunogen information

Immunogen:	FK506 binding protein 9, 63 kDa
Synonyms:	Peptidyl-prolyl cis-trans isomerase FKBP9 (PPIase FKBP9) 63 kDa FK506-binding protein (63 kDa FKBP, FKBP-63) FK506-binding protein 9 (FKBP-9) Rotamase FKBP9 FKBP60 FKBP63
Observed MW:	63 kDa
Uniprot ID :	O95302

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab03150(FKBP9 Antibody) at dilution of 1:600