

FMN1 antibody

Product Information

Catalog No.:	FNab03165
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the formation of adherens junction and the polymerization of linear actin cables.

Immunogen information

Immunogen:	formin 1
Synonyms:	Formin-1 Limb deformity protein homolog FMN1 FMN LD
Observed MW:	160 kDa
Uniprot ID :	Q68DA7

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab03165(FMN1 Antibody) at dilution of 1:600