

## FMO2-specific antibody

### Product Information

Catalog No.:	FNab03169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	flavin containing monooxygenase 2(non-functional)
Synonyms:	Flavin-containing monooxygenase 2 Dimethylaniline oxidase 2 FMO 1B1 Pulmonary flavin-containing monooxygenase 2 (FMO 2) FMO2
Observed MW:	50 kDa
Uniprot ID :	Q99518

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab03169(FMO2-specific antibody) at dilution of 1:500