

## FMR1NB antibody

### Product Information

Catalog No.:	FNab03173
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	fragile X mental retardation 1 neighbor
Synonyms:	FMR1 neighbor protein Cancer/testis antigen 37 (CT37) Sarcoma antigen NY-SAR-35 FMR1NB
Observed MW:	66 kDa
Uniprot ID :	Q8N0W7

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab03173(FMR1NB antibody) at dilution of 1:1000