

FNDC3A antibody

Product Information

Catalog No.:	FNab03176
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Mediates spermatid-Sertoli adhesion during spermatogenesis.

Immunogen information

Immunogen:	fibronectin type III domain containing 3A
Synonyms:	Fibronectin type-III domain-containing protein 3A Human gene expressed in odontoblasts FNDC3A FNDC3 HUGO KIAA0970
Observed MW:	120-130 kDa
Uniprot ID :	Q9Y2H6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03176(FND3A Antibody) at dilution of 1:1000