

## FST antibody

### Product Information


Catalog No.:	FNab03183
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7F10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	follistatin
Synonyms:	Follistatin (FS) Activin-binding protein FST
Observed MW:	35 kDa
Uniprot ID :	P19883

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab03183( Follistatin antibody) at dilution of 1:400



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03183(Follistatin antibody) at dilution of 1:1000