

## FOXC2 antibody

### Product Information

Catalog No.:	FNab03191
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcriptional activator. Might be involved in the formation of special mesenchymal tissues.

### Immunogen information

Immunogen:	forkhead box C2(MFH-1, mesenchyme forkhead 1)
Synonyms:	Forkhead box protein C2 Forkhead-related protein FKHL14 Mesenchyme fork head protein 1 (MFH-1 protein) Transcription factor FKH-14 FOXC2 FKHL14 MFH1
Observed MW:	56 kDa
Uniprot ID :	Q99958

### Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab03191(FOXC2 antibody) at dilution of 1:1000