

## FOXG1 antibody

### Product Information

Catalog No.:	FNab03195
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

Transcription repression factor which plays an important role in the establishment of the regional subdivision of the developing brain and in the development of the telencephalon.


### Immunogen information

Immunogen:	forkhead box G1
Synonyms:	Forkhead box protein G1 Brain factor 1 (BF-1, BF1) Brain factor 2 (BF-2, BF2, hBF-2) Forkhead box protein G1A Forkhead box protein G1B Forkhead box protein G1C Forkhead-related protein FKHL1 (HFK1) Forkhead-related protein FKHL2 (HFK2) Forkhead-related protein FKHL3 (HFK3) FOXG1 FKH2 FKHL1 FKHL2 FKHL3 FKHL4 FOXG1A FOXG1B FOXG1C
Observed MW:	60 kDa
Uniprot ID :	P55316

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Image:



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab03195(FOXG1 antibody) at dilution of 1:500