

## FOXH1-Specific antibody

### Product Information

Catalog No.:	FNab03196
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	forkhead box H1
Synonyms:	Forkhead box protein H1 Forkhead activin signal transducer 1 (Fast-1, hFAST-1) Forkhead activin signal transducer 2 (Fast-2) FOXH1 FAST1 FAST2
Observed MW:	42 kDa
Uniprot ID :	O75593

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:1000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03196(FOXH1-Specific antibody) at dilution of 1:300