

## FOXJ3 antibody

### Product Information

Catalog No.:	FNab03197
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcriptional activator of MEF2C involved in the regulation of adult muscle fiber type identity and skeletal muscle regeneration.

### Immunogen information

Immunogen:	forkhead box J3
Synonyms:	Forkhead box protein J3 FOXJ3 KIAA1041
Observed MW:	66-69 kDa
Uniprot ID :	Q9UPW0

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab03197(FOXJ3 antibody) at dilution of 1:500