

## FOXL2-Specific antibody

### Product Information

Catalog No.:	FNab03199
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	forkhead box L2
Synonyms:	Forkhead box protein L2 FOXL2
Observed MW:	49 kDa
Uniprot ID :	P58012

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab03199( FOXL2-Specific Antibody) at dilution of 1:300