

FOXP1 antibody

Product Information


Catalog No.:	FNab03206
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	forkhead box P1
Synonyms:	Forkhead box protein P1 Mac-1-regulated forkhead (MFH) FOXP1
Observed MW:	50 kDa, 60-65 kDa, 85 kDa
Uniprot ID :	Q9H334

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03206(FOXP1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded Rat spleen using FNab3206(FOXP1 Antibody) at dilution of 1:100