

## FOXQ1 antibody

### Product Information

Catalog No.:	FNab03211
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Plays a role in hair follicle differentiation.

### Immunogen information

Immunogen:	forkhead box Q1
Synonyms:	Forkhead box protein Q1 HNF-3/forkhead-like protein 1 (HFH-1) Hepatocyte nuclear factor 3 forkhead homolog 1 FOXQ1 HFH1
Observed MW:	42 kDa
Uniprot ID :	Q9C009

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab03211(FOXQ1 antibody) at dilution of 1:800