

ATAD5 antibody

Product Information

Catalog No.:	FNab03217
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	FGF receptor activating protein 1
Synonyms:	ATPase family AAA domain-containing protein 5 Chromosome fragility-associated gene 1 protein ATAD5 C17orf41 Elg1 FRAG1
Observed MW:	150 kDa
Uniprot ID :	Q96QE3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03217(ATAD5 Antibody) at dilution of 1:600