

FZD9 antibody

Product Information

Catalog No.:	FNab03223
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	frizzled homolog 9(Drosophila)
Synonyms:	Frizzled-9 (Fz-9, hFz9) FzE6 FZD9 FZD3
Observed MW:	64 kDa
Uniprot ID :	O00144

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab03223(Frizzled 9 antibody) at dilution of 1:2000