

FSD1L antibody

Product Information

Catalog No.:	FNab03229
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	fibronectin type III and SPRY domain containing 1-like
Synonyms:	FSD1-like protein Coiled-coil domain-containing protein 10 FSD1 N-terminal-like protein FSD1L CCDC10 CSDUFD1 FSD1CL FSD1NL
Observed MW:	56-65 kDa
Uniprot ID :	Q9BXM9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab03229(FSD1L antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03229(FSD1L antibody) at dilution of 1:500