

FUT4 antibody

Product Information

Catalog No.:	FNab03248
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

May catalyze alpha-1,3 glycosidic linkages involved in the expression of Lewis X/SSEA-1 and VIM-2 antigens.

Immunogen information

Immunogen:	fucosyltransferase 4(alpha(1,3) fucosyltransferase, myeloid-specific)
Synonyms:	Alpha-(1,3)-fucosyltransferase 4 4-galactosyl-N-acetylglucosaminide 3-alpha-L-fucosyltransferase ELAM-1 ligand fucosyltransferase Fucosyltransferase 4 Fucosyltransferase IV (Fuc-TIV, FucT-IV) Galactoside 3-L-fucosyltransferase FUT4 ELFT FCT3A
Observed MW:	95-110 kDa
Uniprot ID :	P22083

Application

Reactivity:	Human
Tested Application:	ELISA