

TIMM10B antibody

Product Information

Catalog No.:	FNab03254
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	fracture callus 1 homolog(rat)
Synonyms:	Mitochondrial import inner membrane translocase subunit Tim10 B Fracture callus protein 1 FXC1 Mitochondrial import inner membrane translocase subunit Tim9 B TIMM10B (Tim10b) TIMM10B FXC1 TIM9B TIMM9B
Observed MW:	None
Uniprot ID :	Q9Y5J6

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA