

FXR2 antibody

Product Information

Catalog No.:	FNab03257
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


RNA-binding protein.

Immunogen information


Immunogen:	fragile X mental retardation, autosomal homolog 2
Synonyms:	RNA-binding protein FXR2 (FXR2P) FMR1 autosomal homolog 2 FXR2 FMR1L2
Observed MW:	95 kDa
Uniprot ID :	P51116

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-FXR2 (IP:FNab03257, 4ug;
Detection:FNab03257 1:500) with HeLa cells
lysate 3200ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab03257(FXR2 antibody)
at dilution of 1:500