

## G0S2 antibody

### Product Information

Catalog No.:	FNab03266
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.

### Immunogen information

Immunogen:	G0/G1switch 2
Synonyms:	G0/G1 switch protein 2 G0/G1 switch regulatory protein 2 Putative lymphocyte G0/G1 switch gene G0S2
Observed MW:	11 kDa
Uniprot ID :	P27469

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab03266(G0S2 antibody) at dilution of 1:1000