

## GA repeat antibody

### Product Information

Catalog No.:	FNab03270
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	None
Synonyms:	None
Observed MW:	None
Uniprot ID :	None

### Application

Reactivity:	Human
Tested Application:	ELISA