

GAB3-Specific antibody

Product Information

Catalog No.:	FNab03273
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	GRB2-associated binding protein 3
Synonyms:	GRB2-associated-binding protein 3 GRB2-associated binder 3 Growth factor receptor bound protein 2-associated protein 3 GAB3
Observed MW:	63 kDa
Uniprot ID :	Q8WWW8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03273(GAB3-Specific Antibody) at dilution of 1:200



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03273(GAB3-Specific antibody) at dilution of 1:500