

GABRA4 antibody

Product Information

Catalog No.:	FNab03285
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Immunogen information

Immunogen:	gamma-aminobutyric acid(GABA) A receptor, alpha 4
Synonyms:	Gamma-aminobutyric acid receptor subunit alpha-4 GABA(A) receptor subunit alpha-4 GABRA4
Observed MW:	62-70 kDa
Uniprot ID :	P48169

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab03285(GABRA4 antibody) at dilution of 1:500