

LGALS4 antibody

Product Information

Catalog No.:	FNab03311
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	lectin, galactoside-binding, soluble, 4
Synonyms:	Galectin-4 (Gal-4) Antigen NY-CO-27 L-36 lactose-binding protein (L36LBP) Lactose-binding lectin 4 LGALS4
Observed MW:	38 kDa
Uniprot ID :	P56470

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab03311(LGALS4,GAL4 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03311(LGALS4,GAL4 antibody) at dilution of 1:400