

LGALS8 antibody

Product Information

Catalog No.:	FNab03312
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	lectin, galactoside-binding, soluble, 8
Synonyms:	Galectin-8 (Gal-8) Po66 carbohydrate-binding protein (Po66-CBP) Prostate carcinoma tumor antigen 1 (PCTA-1) LGALS8
Observed MW:	None
Uniprot ID :	O00214

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA