

LGALS1 antibody

Product Information


Catalog No.:	FNab03315
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9F4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	lectin, galactoside-binding, soluble, 1
Synonyms:	Galectin-1 (Gal-1) 14 kDa laminin-binding protein (HLBP14) 14 kDa lectin Beta-galactoside-binding lectin L-14-I Galaptin HBL HPL Lactose-binding lectin 1 Lectin galactoside-binding soluble 1 Putative MAPK-activating protein PM12 S-Lac lectin 1 LGALS1
Observed MW:	15 kDa
Uniprot ID :	P09382

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab03315(GAL1 antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03315(GAL1 antibody) at dilution of 1:1000