

## LGALS7 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab03318  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information


|              |  |
|--------------|--|
| Immunogen:   | lectin, galactoside-binding, soluble, 7B                             |
| Synonyms:    | Galectin-7 (Gal-7) HKL-14 PI7 p53-induced gene 1 protein LGALS7 PIG1 |
| Observed MW: | 15 kDa   |
| Uniprot ID : | P47929   |

### Application

|                       |                                    |
|-----------------------|------------------------------------|
| Reactivity:           | Human, Mouse, Rat                  |
| Tested Application:   | ELISA, WB, IHC                     |
| Recommended dilution: | WB: 1:1000-1:4000; IHC: 1:20-1:200 |
| Image:                |                                    |



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab03318(Galectin-7 antibody) at dilution of 1:100



mouse skin tissue were subjected to SDS PAGE followed by western blot with FNab03318(Galectin-7 antibody) at dilution of 1:1500