

GALNT15 antibody

Product Information

Catalog No.:	FNab03327
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None


Immunogen information

Immunogen:	UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase-like 2
Synonyms:	Polypeptide N-acetylgalactosaminyltransferase 15 Polypeptide GalNAc transferase-like protein 2 (GalNAc-T-like protein 2, pp-GaNTase-like protein 2) Polypeptide N-acetylgalactosaminyltransferase-like protein 2 Protein-UDP acetylgalactosaminyltransferase-like protein 2 UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase-like protein 2 GALNT15 GALNTL2
Observed MW:	65-70 kDa
Uniprot ID :	Q8N3T1


Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human placenta using FNab03327(GALNTL2 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03327(GALNTL2 antibody) at dilution of 1:200