

AP1G1 antibody

Product Information

Catalog No.:	FNab03332
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	adaptor-related protein complex 1, gamma 1 subunit
Synonyms:	AP-1 complex subunit gamma-1 Adaptor protein complex AP-1 subunit gamma-1 Adaptor-related protein complex 1 subunit gamma-1 Clathrin assembly protein complex 1 gamma-1 large chain Gamma1-adaptin Golgi adaptor HA1/AP1 adaptin subunit gamma-1 AP1G1 ADTG CLAPG1
Observed MW:	91 kDa
Uniprot ID :	O43747

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab03332 (Gamma adaptin Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03332 (Gamma adaptin antibody) at dilution of 1:400