

YEATS4 antibody

Product Information

Catalog No.:	FNab03351
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	YEATS domain containing 4
Synonyms:	YEATS domain-containing protein 4 Glioma-amplified sequence 41 (Gas41) NuMA-binding protein 1 (NuBI-1, NuBI1) YEATS4 GAS41
Observed MW:	27-30 kDa
Uniprot ID :	O95619

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03351(GAS41 antibody) at dilution of 1:300