

GAS6 antibody

Product Information

Catalog No.:	FNab03352
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a gamma-carboxyglutamic acid (Gla)-containing protein thought to be involved in the stimulation of cell proliferation. This gene is frequently overexpressed in many cancers and has been implicated as an adverse prognostic marker. Elevated protein levels are additionally associated with a variety of disease states, including venous thromboembolic disease, systemic lupus erythematosus, chronic renal failure, and preeclampsia.

Immunogen information


Immunogen:	growth arrest-specific 6
Synonyms:	Growth arrest-specific protein 6 (GAS-6) AXL receptor tyrosine kinase ligand GAS6 AXLLG
Observed MW:	80 kDa
Uniprot ID :	Q14393

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse lung using FNab03352(GAS6 antibody) at dilution of 1:50



DU 145 cells were subjected to SDS PAGE followed by western blot with FNab03352(GAS6 Antibody) at dilution of 1:1000