

GAS7 antibody

Product Information

Catalog No.:	FNab03353
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in promoting maturation and morphological differentiation of cerebellar neurons.

Immunogen information

Immunogen:	growth arrest-specific 7
Synonyms:	Growth arrest-specific protein 7 (GAS-7) GAS7 KIAA0394
Observed MW:	38-50 kDa
Uniprot ID :	O60861

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab03353(GAS7 Antibody) at dilution of 1:1000