

GAST antibody

Product Information

Catalog No.:	FNab03355
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	gastrin
Synonyms:	Gastrin Gastrin-71 Alternative names: Gastrin component I Gastrin-52 (G52) Big gastrin Alternative names: Gastrin component II Gastrin-34 (G34) Gastrin Alternative names: Gastrin component III Gastrin-17 (G17) Gastrin-14 (G14) Gastrin-6 (G6) GAST GAS
Observed MW:	11 kDa
Uniprot ID :	P01350

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab03355(Gastrin antibody) at dilution of 1:200