

GKN1 antibody

Product Information

Catalog No.:	FNab03356
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	gastrokine 1
Synonyms:	Gastrokine-1 18 kDa antrum mucosa protein (AMP-18) Protein CA11 GKN1 AMP18 CA11
Observed MW:	20-22 kDa
Uniprot ID :	Q9NS71

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab03356(GKN1 antibody) at dilution of 1:100



human stomach tissue were subjected to SDS PAGE followed by western blot with FNab03356(GKN1 antibody) at dilution of 1:2000