

## GATA4 antibody

### Product Information

Catalog No.:	FNab03363
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	GATA binding protein 4
Synonyms:	Transcription factor GATA-4 GATA-binding factor 4 GATA4
Observed MW:	45 kDa
Uniprot ID :	P43694

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab03363(GATA4-Specific antibody) at dilution of 1:50



IP Result of anti-GATA4-Specific (IP:FNab03363, 4ug; Detection:FNab03363 1:600) with mouse heart tissue lysate 3200ug.



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03363(GATA4-Specific antibody) at dilution of 1:500