

## GATA5 antibody

### Product Information

Catalog No.:	FNab03364
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor required during cardiovascular development(PubMed:23289003). Plays an important role in the transcriptional program(s) that underlies smooth muscle cell diversity(By similarity). Binds to the functionally important CEF-1 nuclear protein binding site in the cardiac-specific slow/cardiac troponin C transcriptional enhancer(PubMed:25543888).

### Immunogen information


Immunogen:	GATA binding protein 5
Synonyms:	Transcription factor GATA-5 GATA-binding factor 5 GATA5
Observed MW:	41 kDa
Uniprot ID :	Q9BWX5

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab03364(GATA5 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03364(GATA5 Antibody) at dilution of 1:300