

GBP6 antibody

Product Information

Catalog No.:	FNab03377
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds GTP, GDP and GMP.

Immunogen information


Immunogen:	guanylate binding protein family, member 6
Synonyms:	Guanylate-binding protein 6 GTP-binding protein 6 (GBP-6) Guanine nucleotide-binding protein 6 GBP6
Observed MW:	70 kDa
Uniprot ID :	Q6ZN66

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab03377(GBP6 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03377(GBP6 Antibody) at dilution of 1:300