

GUCA1A antibody

Product Information

Catalog No.:	FNab03378
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	guanylate cyclase activator 1A(retina)
Synonyms:	Guanylyl cyclase-activating protein 1 (GCAP 1) Guanylate cyclase activator 1A GUCA1A C6orf131 GCAP GCAP1 GUCA1
Observed MW:	23-24 kDa
Uniprot ID :	P43080

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Y79 cells were subjected to SDS PAGE followed by western blot with FNab03378(GCAP1 antibody) at dilution of 1:200