

## GCC1 antibody

### Product Information

Catalog No.:	FNab03380
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Probably involved in maintaining Golgi structure.

### Immunogen information

Immunogen:	GRIP and coiled-coil domain containing 1
Synonyms:	GRIP and coiled-coil domain-containing protein 1 Golgi coiled-coil protein 1 GCC1
Observed MW:	105 kDa
Uniprot ID :	Q96CN9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03380( GCC1 Antibody) at dilution of 1:2000